Advancis Surgical

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CE 1639

Advancis Surgical is a trading name of Brightwake Limited.

Advancis Surgical

Hemosep® Cell Concentrator Pack

- Directions for use**
1. Open Hemosep Cell Concentrator Pack.
 2. Remove cell concentrator bag from packaging.
 3. Close the male and one of the female ports using the clips.
 4. Prime the cell concentrator bag by injecting 100ml saline (0.9% w/v) into the open female port using a luer lock syringe or giving set.
 5. Close the clip.
 6. Rock the cell concentrator bag for 30 seconds. Do not squeeze or bend bag/ or damage the membrane.
 7. The membrane should be a light grey colour. Blood must be introduced within 3 minutes of priming.
 8. Connect Hemosep Cell Concentrator Bag to Heart Lung Machine Reservoir or collection reservoir via either of the inlet ports.
 9. Open the clip & fill the bag to a max 500ml per bag (>500ml will decrease PCV end result).
 10. Close the clip and place on the Hemosep Shaker.
 11. Press run and process for 20 minutes. Reducing this time will decrease the PCV end result.
 12. Remove the Hemosep Cell Concentrator Bag from the Hemosep Shaker.
 13. Connect Hemosep Cell Concentrator Bag to Hemosep Blood Collection Bag.
 14. Open the clip and hold the Hemosep Cell Concentrator Bag above the Hemosep Blood Collection Bag.
 15. Allow the processed blood to drain into the Hemosep Blood Collection Bag by gravity.
 16. Close clip.
 17. Disconnect the Hemosep Cell Concentrator Bag from the Hemosep Blood Collection Bag.
 18. The blood is now ready for transfusion through a standard blood giving set with a 170-260 micron filter.
 19. Dispose of Hemosep Cell Concentrator Bag in clinical waste.

Indications for use

The use of Hemosep for intraoperative cell salvage is recommended in a variety of medium to high blood loss procedures including cardiac, orthopaedic, spinal, renal and vascular surgeries. Hemosep cell salvage may be used in surgical procedures where the surgical field is not contaminated by faecal or infective matter and where no other contraindications exist. Patient selection for Hemosep cell salvage is at the discretion of the surgeon and anaesthetist caring for the patient. Each organisation should refer to its local policy for use of intraoperative cell salvage.

Contraindications

Hemosep cell salvage should not be used in the following situations:

- Bowel contents in the surgical field (the surgical field may be irrigated with copious normal saline (0.9% NaCl) before resuming Hemosep cell salvage).
- During obstetric surgery due to the risk of amniotic fluid in the surgical field.
- Heparin induced thrombocytopenia if heparin is the only available anticoagulant (a citrate anticoagulant solution may be used instead).

Cautions

- Pleural effusions should not be aspirated and should be drained prior to Hemosep cell salvage.
- The decision to use Hemosep cell salvage in patients with abnormal red cell disorders, such as sickle cell disease, should be made by a clinician on an individual patient basis.
- The use of Hartmann's Solution will inhibit the action of citrate based anticoagulants (e.g. ACD) if used as an irrigant.
- Where possible, all air in the Hemosep blood collection should be evacuated prior to reinfusion. Do not use a pressure cuff as there is a risk of air embolus.
- The decision to use blood that is potentially contaminated with malignant cells should be made by the clinicians caring for the patient, taking into account the latest evidence, considering the risk and benefits for the individual patient. Clinicians should consider using a leukoreduction filter to reinfuse the blood.¹
- During high blood loss periods, it may become necessary to increase the suction at the request of the surgical team. The vacuum should be returned to standard level (100-150mmHg) as soon as bleeding is under control. This will minimise damage to the RBC and haemolysis.¹
- Where possible the tip of the wand should be immersed in blood and NOT skimmed across the surface of tissues or pools of blood. Skimming results in a large quantity of air mixing with the aspirated blood, which increases haemolysis and therefore reduces the number of viable RBC for reinfusion.¹
- When substances not licensed for intravenous (IV) use are used within the surgical field and could potentially be aspirated into the collection reservoir: In this instance, a separate standard theatre suction should be used to remove the substance. The surgical field should then be irrigated with copious normal saline (0.9% NaCl) before resuming Hemosep cell salvage.
- Examples of non-IV materials that should not be aspirated into the Hemosep cell salvage system include: antibiotics not licensed for IV use, iodine, topical clotting agents, methyl methacrylate (orthopaedic bone cement).
- The presence of infection may result in contamination of the salvaged blood. The decision to use blood that is potentially contaminated with bacteria should be made by the clinicians caring for the patient, taking into account the latest evidence, considering the risk and benefits for the individual patient.
- Where gastric/pancreatic secretions could be aspirated into the system as they may cause enzymatic haemolysis.

References

1. Section 8 – Intraoperative Cell Salvage Education Workbook. <https://www.transfusionguidelines.org/transfusion-practice/uk-cell-salvage-action-group/intraoperative-cell-salvage-education>

Dispose

Dispose as 'Clinical Waste'.

Device Accessory List

The Hemosep® system comprises of the following parts:

- Hemosep® Cell Concentrator Pack, Hemosep® Shaker Unit, Hemosep® Intra-Operative stand, Hemosep® Intra-Operative Kit.

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