


Eclipse NON-BACKED

The maximum absorbency superabsorbent dressing

- 1 Wound Contact layer
- 2 Double CrystaLock layer
- 3 Rapid wicking layer
- 4 Soft edges for comfort





Eclipse Non-Backed absorbs an extremely high level of exudate

- Reduced need for stacking
- CrystaLock technology locks away exudate, bacteria and MMPs²
- Central wicking layer for even distribution of exudate



Eclipse NON-BACKED

Eclipse Non-Backed offers market leading absorbency



Ordering Information

Size	Product Code	Pack Size	PIP Code
10 x 10cm	CR4516	20	419-7836
10 x 20cm	CR4517	10	419-7844
20 x 20cm	CR4518	10	419-7935
20 x 30cm	CR4519	20	419-7927

Get in touch

+44 (0)1623 751 500
info@advancis.co.uk

Find us online

www.advancis.co.uk

Advancis, Eclipse and CrystalLock are registered trademarks of Brightwake Ltd. Contact info@advancis.co.uk for more information

References:

1. Rafter, L. Anthony, D. Collier, M. Rafter, R. (2015) Stopping the strikethrough: An audit of patient outcomes on four superabsorbent dressings. Wounds UK Vol 11 No 4, 60-67
2. Weigand, C. White, R. (2013) Binding and inhibition of protease enzymes, including MMPs, by a superabsorbent in vitro. Journal of Wound Care Vol 22, No 5