

# Siltape®

Soft silicone perforated tape

## Preventing facial skin damage beneath PPE with Siltape

During the **COVID-19 outbreak**, healthcare professionals are seeing early signs of **skin damage** with the increased use of facial masks. Accurate fit testing and the use of barrier creams can help. If additional protection is needed under a mask, **Siltape can be used.\***

\*When using fit tested masks, ensure the fit test is done after the application of Siltape and that the seal isn't compromised.

## What is Siltape?

Siltape is a soft silicone tape on a handy roll with perforations for ease of use; simply tear off what you need and apply to the affected area.

## Where to use Siltape:

- Bridge of nose
- Cheeks
- Tops of ears
- Chin
- Forehead



## Ordering Information

Size	Stock code	Pack size	NHS code	PIP code
2cm x 3m	CR3938	1	EHR053	344-0336
4cm x 1.5m	CR3939	1	EHR054	344-0328